

Technical Data



Seat width
34 – 44 cm
(in increments
of 2 cm)



Seat depth
37.5 – 45 cm
(in increments
of 2.5 cm)



Seat height
front: 49 – 51 cm
rear: 40 – 49 cm
(in increments
of 1 cm)



Backrest height:
27 – 48 cm
(in increments
of 1.5 cm)



Knee-to-heel length:
39 – 48 cm
(in increments
of 1 cm)



Seat angle:
variable,
standard 6 cm



Backrest angle:
-12°, -8°, -4°, 0°, 4°,
8°, 12°



Total width with
3° camber:
Seat width + 22.5 cm



Total length:
75°: approx. 83 cm;
90°: approx. 76 cm



Total weight/
unloaded:
approx. 9.8 kg (with
seat width 36 cm)



Max. user weight:
100 kg



Transport weight
(w/out rear wheels):
approx. 6.6 kg (with
seat width 36 cm)

NOTE: Unless otherwise specified, all measurements and weights indicated are based on the standard configuration of this model.

Frame colours



Glossy black



Night shimmer



Silver grey



Mirror natural



Royal red

CPS (Colour Power Set)



Techno anthracite



Lemon



Comet red



Arctic blue



B-blue

CPS includes the "Starec" hubs, 24" rims, back angle adjustment plate, centre part, forks as well as castor housing and seat module caps.

Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Tel +41 (0)61 487 70 70 | Fax +41 (0)61 487 70 79 | www.kuschallR33.com

United Kingdom:
Invacare Ltd. | South Road | Bridgend Industrial Estate | Bridgend CF31 3PY | Phone +44 (0)1 656 647327 | Fax +44 (0)1 656 649016

Ireland:
Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085

Germany and Eastern Europe:
Invacare Deutschland GmbH | Kleiststrasse 49 | D-32457 Porta Westfalica | Tel +49 (0) 5731 7540 | Fax +49 (0) 5731 754 52191

Küschall® is a registered trademark.
Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication, or of parts thereof, as well as any reproduction of images, must have the written consent of Küschall AG.



Designer technology for the resolute

My dreams are to make yours reality. With the *kuschall R33*, I combined state-of-the-art technology and revolutionary design to create a wheelchair for the champion within all of us. The *R33* is passion: my passion for success and victory.

R. Kuschall
Rainer Kuschall

Inspired Technology

kuschall **R33**



The future starts with today's technology

Experience the new revolutionary shock absorbing system with this state-of-the-art wheelchair.

More than just a wheelchair

Because style is just as important as performance.



The difference lies in detail

A powerful design for performance leaves a lasting impression.